

THIS PACKAGE INSERT CONTAINS IMPORTANT PRODUCT USE AND SAFETY INFORMATION. PLEASE READ CAREFULLY AND RETAIN FOR FUTURE REFERENCE.



**One Bottle Solution
for Cleaning & Disinfecting**

DESCRIPTION / CONTENTS:

Clear Care Cleaning & Disinfecting Solution is a sterile solution containing micro-filtered hydrogen peroxide 3%, sodium chloride 0.79%, stabilized with phosphonic acid, a phosphate buffered system, and Pluronic 17R4 (a cleaning agent).

GOOD LENS CARE PRACTICES :

To ensure proper disinfection of your lenses you must follow the instructions completely. Do not skip any steps. Always wash and rinse your hands before handling your lenses. This will help prevent eye infections by removing dirt and oils that could get on the lenses.

DIRECTIONS FOR USE - Soft Lenses:


Do not put Clear Care Cleaning and Disinfecting Solution on your lenses and insert directly into the eye.

TO CLEAN, DISINFECT, NEUTRALIZE, AND REMOVE PROTEIN FROM YOUR LENSES

- Remove and place each lens into the appropriately marked L/R domed lens holder. Rinse with Clear Care for 5 seconds.
- Fill the lens case to fill line with Clear Care and place the lens holder in the case.
- Tighten the cap and store lenses for at least 6 hours or overnight. **DO NOT SHAKE THE CASE.**

NOTE: To prevent damage to your lens, center the lens on the dome in the lens holder.

Be sure the lens does not touch the basket rim, then close the basket lid.



After soaking for six hours, your lenses are ready to wear. No final saline rinse is necessary. **Never rinse your lenses with Clear Care Cleaning and Disinfecting Solution prior to insertion.** If desired, lenses can be rinsed with a sterile saline before inserting. Discard the neutralized disinfectant from the cup. Rinse the lens cup with fresh saline or Clear Care and allow the case to air dry with the lens holder inverted outside the case. Do not place the lens holder on its side.

Your eye care professional may recommend additional products (such as Unizyme® Enzymatic Cleaner) or procedures to care for your lenses (such as rubbing with PURE EYES® Cleaner/Rinse Solution) based on individual tear chemistry and lens wearing schedule. Always follow your eye care professional's instructions.

- If your eye care professional recommends use of Unizyme Enzymatic Cleaner, follow the directions for use on the Unizyme Enzymatic Cleaner package insert.
- If your eye care professional recommends an additional rubbing step, place two drops of Pure Eyes® Cleaner Rinse on your lens and rub for 20 seconds, then follow the directions for use above.

DIRECTIONS FOR USE - RGP Lenses:

Do not put Clear Care Cleaning and Disinfecting Solution on your lenses and insert directly into the eye.

- **Remove your lenses one at a time and place them into the appropriately marked dome basket holder.**
- Place each lens in the palm of your hand, apply 2 to 4 drops of Clear Care and rub. **While rubbing your lenses with Clear Care, some users may experience a mild, temporary skin discoloration (bleaching) of the fingers or hands. Always wash and rinse your hands after rubbing your lenses with the solution.**
- Return the lenses to the appropriate holder and close the baskets. Thoroughly rinse the lenses for 5 seconds through the basket with Clear Care.
- Fill the lens case with Clear Care and place the lens holder in the case. Tighten the cap and store lenses for at least six (6) hours. **DO NOT SHAKE THE CASE.** Do not rinse the lenses. Place the lenses directly on the eye from the solution or place a few drops of Focus Lens Drops or Aquify Long-Lasting Comfort Drops on the lens for extra cushioning.

STORAGE:

If you do not intend to wear your lenses immediately after disinfection/neutralization, you may store them in the unopened lens case for up to 7 days. Disinfect and neutralize once a week with fresh solution before wearing your lenses.

ACTIONS:

When used as directed Clear Care provides a unique cleaning action, which removes film and debris from the lens surface. Clear Care also helps prevent serious eye infections by killing harmful microorganisms on contact lenses.

INDICATIONS (Uses):

For simultaneous cleaning, daily protein removal, disinfecting, and storing of soft (hydrophilic) or rigid gas permeable (fluoro silicone acrylate and silicone acrylate) contact lenses as recommended by your eye care practitioner.

CONTRAINDICATIONS (Reasons not to use):

There are no known contraindications for use of Clear Care Cleaning & Disinfecting Solution; however, if you are allergic to any ingredient in this solution, do not use.

WARNINGS:

Do not put Clear Care Cleaning and Disinfecting Solution on your lenses and insert directly into the eye. Lenses must be soaked in the provided lens case for 6 hours (neutralization process) prior to lens insertion. Otherwise, burning and stinging will result.

While rubbing RGP lenses with Clear Care, some patients may experience a mild, temporary skin discoloration (bleaching) of the fingers or hands. Always wash and rinse your hands after rubbing your lenses with the solution.

DO NOT PUT CLEAR CARE SOLUTION THAT HAS NOT BEEN NEUTRALIZED IN YOUR EYE. Should unneutralized Clear Care mistakenly get in your eye, remove your lenses immediately, flush (wash) your eyes with a large amount of water or sterile saline for a few minutes. If burning and/or irritation persist, seek assistance from an eye care professional. The red dropper tip indicates that Clear Care Cleaning & Disinfecting Solution should not be put directly in your eye.

PROBLEMS WITH CONTACT LENSES AND LENS CARE PRODUCTS COULD RESULT IN SERIOUS INJURY TO THE EYE.

Follow your eye care professional's directions and all labeling instructions for proper use and care of your lenses and lens care products, including the lens case. Eye problems, including corneal ulcers and infections, can develop rapidly and lead to loss of vision. Daily wear lenses are not indicated for overnight wear and should not be worn while sleeping. Clinical studies have shown the risk of serious adverse reactions is increased when these lenses are worn overnight. Extended wear lenses should be regularly removed for cleaning and disinfection or for disposal and replacement on the schedule prescribed by your eye care professional. Clinical studies have shown that there is an increased incidence of serious adverse reactions in extended wear contact lens users as compared to daily wear contact lens users. Studies have also shown that the risk of serious adverse reactions increases the longer extended wear lenses are worn before removal for cleaning and disinfection or for disposal and replacement. Studies have also shown that smokers have a higher incidence of adverse reactions. If you experience eye discomfort, excessive tearing, vision changes, or redness of the eye, immediately remove your lenses and promptly contact your eye care practitioner. All contact lens wearers should see their eye care practitioner as directed. Periodic eye exams are extremely important in order to detect and treat problems related to contact lens wear before they cause discomfort or pain.

Replace cap after using. To avoid contaminating your solution, do not transfer to other bottles or containers.

All contact lens wearers must see their eye care professional regularly and as frequently as directed. If your lenses are for extended wear, your eye care professional may prescribe more frequent visits.

Studies have shown that smoking increases the risk of ulcerative keratitis for contact lens users.

Do not take Clear Care Cleaning and Disinfecting Solution internally or gastric distress will result. Seek professional assistance of physician or a poison control center immediately if taken internally. Keep out of the reach of children.

PRECAUTIONS:

- **Only use the Clear Care Lens Case for disinfection. Do not use a flat case.**
- Always wash and dry your hands before handling your lenses.
- Never use this for heat disinfection.
- Never reuse this solution.
- Keep bottle tightly closed when not in use.
- Store at room temperature. (15°C to 30°C / 59°F to 86°F)
- Use before the expiration date marked on the container.

There are no safe, acceptable substitutes for Clear Care Cleaning & Disinfecting Solution components. Do not mix or substitute other lens care hydrogen peroxide products or lens cases. Do not use any other disinfection solutions in the disposable lens case. **DO NOT USE OVER-THE-COUNTER GENERIC HYDROGEN PEROXIDE.** Generic hydrogen peroxide solutions are not intended for ophthalmic use and may contain ingredients not tested for ocular safety or toxicity. Generic hydrogen peroxide may contain ingredients that cause **DISCOLORATION OR DAMAGE YOUR CONTACT LENSES.**

To avoid possible damage to your lenses, the following precautions should be taken when using Clear Care Cleaning & Disinfecting Solution:

- Do not over tighten the case (only tighten finger tight).
- Do not use a damaged lens case.
- Do not fill above the fill line on the case.
- Do not shake the case during disinfection.
- Dispose of your old lens case with each new purchase of Clear Care Cleaning and Disinfecting Solution. A fresh new lens case is supplied with each purchase.

ADVERSE REACTIONS (Possible problems and what to do):

The following problems may occur with contact lens wear: eyes sting, burn or itch (irritation), comfort is less than when lens was first placed on the eye, feeling of something in the eye (foreign body, scratched area), excessive watering (tearing) of the eye, unusual eye secretions, redness of the eye, reduced sharpness of vision (poor visual acuity), blurred vision, rainbows or halos around objects, sensitivity to light (photophobia), or dry eyes.

If you notice any of the above, **IMMEDIATELY** remove your lenses.

- If the problem stops and the lenses appear to be undamaged, repeat the cleaning, disinfection and neutralization process and reinsert them. **Do not put unneutralized solution in your eyes.**
- If the lens is in any way damaged, do not put the lens back on your eye. Place the lens in the storage case and contact your eye care practitioner.
- If the lens has dirt, an eyelash, or other foreign body on it, repeat the cleaning, disinfection and neutralization process and reinsert them. **DO NOT PUT Clear Care Solution that has not been neutralized in your eye.**
- If the problem continues **IMMEDIATELY** remove the lens and consult your eye care professional.

If any of the above symptoms occur, a serious condition such as infection, corneal ulcer, neovascularization or iritis may be present. Seek immediate professional identification of the problem and prompt treatment to avoid serious eye damage.

Eye Care Professionals: All adverse reactions observed using the Clear Care Cleaning and Disinfecting Solution should be reported to:

Consumer Relations, CIBA Vision Corporation, Duluth, GA 30097. 1-800-875-3001

HOW SUPPLIED:

- Clear Care Cleaning and Disinfecting Solution is supplied in sterile 4 fl. oz. (118 ml), 8 fl. oz. (237 ml) and 12 fl. oz. (355 ml) plastic bottles. These bottles are marked with lot number and expiration date.
- **TAMPER RESISTANT:** Do not use if CIBA Vision Safety Seal around the bottle cap is broken or missing.

Protected by U.S. Patents 4,889,689; 5,558,846; 5,609,264; 5,609,837 and other patents pending.

For more information regarding this product consult your eye care professional, or write to:

CIBA Vision Corporation, Consumer Relations, Duluth, GA 30097

©2004 CIBA Vision Corporation

Made in Canada

Manufactured for:

CIBA Vision Corporation, Duluth, GA 30097 U.S.A.

3/2005
S7422A



