

Focus® Progressives

Initial Lens Power Selection Tables

$$\text{Initial Power} = \text{Spherical equivalent refraction} + \frac{\text{spherical Add}}{2}$$

(vertex corrected)

eg. SRx = -2.00; +1.50 Add: Trial lens power = -1.25D


MYOPIA Equivalent Sphere Spectacle Refraction*

| | ADD | | | |
|--------|--------------|-------|--------------|--------------|
| | 1.00 to 1.25 | 1.50 | 1.75 to 2.25 | 2.50 to 3.00 |
| 0.00 | +0.50 | +0.75 | +1.00 | +1.25 |
| -0.25 | +0.25 | +0.50 | +0.75 | +1.00 |
| -0.50 | 0.00 | +0.25 | +0.50 | +0.75 |
| -0.75 | -0.25 | 0.00 | +0.25 | +0.50 |
| -1.00 | -0.50 | -0.25 | 0.00 | +0.25 |
| -1.25 | -0.75 | -0.50 | -0.25 | 0.00 |
| -1.50 | -1.00 | -0.75 | -0.50 | -0.25 |
| -1.75 | -1.25 | -1.00 | -0.75 | -0.50 |
| -2.00 | -1.50 | -1.25 | -1.00 | -0.75 |
| -2.25 | -1.75 | -1.50 | -1.25 | -1.00 |
| -2.50 | -2.00 | -1.75 | -1.50 | -1.25 |
| -2.75 | -2.25 | -2.00 | -1.75 | -1.50 |
| -3.00 | -2.50 | -2.25 | -2.00 | -1.75 |
| -3.25 | -2.75 | -2.50 | -2.25 | -2.00 |
| -3.50 | -3.00 | -2.75 | -2.50 | -2.25 |
| -3.75 | -3.25 | -3.00 | -2.75 | -2.50 |
| -4.00 | -3.50 | -3.25 | -3.00 | -2.75 |
| -4.25 | -3.50 | -3.25 | -3.00 | -2.75 |
| -4.50 | -3.75 | -3.50 | -3.25 | -3.00 |
| -4.75 | -4.00 | -3.75 | -3.50 | -3.25 |
| -5.00 | -4.25 | -4.00 | -3.75 | -3.50 |
| -5.25 | -4.50 | -4.25 | -4.00 | -3.75 |
| -5.50 | -4.75 | -4.50 | -4.25 | -4.00 |
| -5.75 | -5.00 | -4.75 | -4.50 | -4.25 |
| -6.00 | -5.25 | -5.00 | -4.75 | -4.50 |
| -6.25 | -5.25 | -5.00 | -4.75 | -4.50 |
| -6.50 | -5.50 | -5.25 | -5.00 | -4.75 |
| -6.75 | -5.75 | -5.50 | -5.25 | -5.00 |
| -7.00 | -6.00 | -5.75 | -5.50 | -5.25 |
| -7.25 | -6.25 | -6.00 | -5.75 | -5.50 |
| -7.50 | -6.50 | -6.25 | -6.00 | -5.75 |
| -7.75 | -6.75 | -6.50 | -6.25 | -6.00 |
| -8.00 | -7.00 | -6.75 | -6.50 | -6.25 |
| -8.25 | -7.25 | -7.00 | -6.75 | -6.50 |
| -8.50 | -7.25 | -7.00 | -6.75 | -6.50 |
| -8.75 | -7.50 | -7.25 | -7.00 | -6.75 |
| -9.00 | -7.75 | -7.50 | -7.25 | -7.00 |
| -9.25 | -8.00 | -7.75 | -7.50 | -7.25 |
| -9.50 | -8.25 | -8.00 | -7.75 | -7.50 |
| -9.75 | -8.25 | -8.00 | -7.75 | -7.50 |
| -10.00 | -8.50 | -8.25 | -8.00 | -7.75 |
| -10.25 | -8.75 | -8.50 | -8.25 | -8.00 |
| -10.50 | -8.75 | -8.50 | -8.25 | -8.00 |
| -10.75 | -9.00 | -8.75 | -8.50 | -8.25 |
| -11.00 | -9.00 | -9.00 | -8.75 | -8.50 |

HYPEROPIA Equivalent Sphere Spectacle Refraction*

| | ADD | | | |
|-------|--------------|-------|--------------|--------------|
| | 1.00 to 1.25 | 1.50 | 1.75 to 2.25 | 2.50 to 3.00 |
| 0.00 | +0.50 | +0.75 | +1.00 | +1.25 |
| +0.25 | +0.75 | +1.00 | +1.25 | +1.50 |
| +0.50 | +1.00 | +1.25 | +1.50 | +1.75 |
| +0.75 | +1.25 | +1.50 | +1.75 | +2.00 |
| +1.00 | +1.50 | +1.75 | +2.00 | +2.25 |
| +1.25 | +1.75 | +2.00 | +2.25 | +2.50 |
| +1.50 | +2.00 | +2.25 | +2.50 | +2.75 |
| +1.75 | +2.25 | +2.50 | +2.75 | +3.00 |
| +2.00 | +2.50 | +2.75 | +3.00 | +3.25 |
| +2.25 | +2.75 | +3.00 | +3.25 | +3.50 |
| +2.50 | +3.00 | +3.25 | +3.50 | +3.75 |
| +2.75 | +3.25 | +3.50 | +3.75 | +4.00 |
| +3.00 | +3.50 | +3.75 | +4.00 | +4.25 |
| +3.25 | +3.75 | +4.00 | +4.25 | +4.50 |
| +3.50 | +4.00 | +4.25 | +4.50 | +4.75 |
| +3.75 | +4.25 | +4.50 | +4.75 | +5.00 |
| +4.00 | +4.75 | +5.00 | +5.25 | +5.50 |
| +4.25 | +5.00 | +5.25 | +5.50 | +5.75 |
| +4.50 | +5.25 | +5.50 | +5.75 | +6.00 |
| +4.75 | +5.50 | +5.75 | +6.00 | +6.25 |
| +5.00 | +5.75 | +6.00 | +6.25 | +6.50 |
| +5.25 | +6.00 | +6.25 | +6.50 | +6.75 |
| +5.50 | +6.25 | +6.50 | +6.75 | +7.00 |
| +5.75 | +6.50 | +6.75 | +7.00 | +7.25 |
| +6.00 | +7.00 | +7.25 | +7.50 | +7.75 |
| +6.25 | +7.25 | +7.50 | +7.75 | +8.00 |
| +6.50 | +7.50 | +7.75 | +8.00 | +8.00 |
| +6.75 | +7.75 | +8.00 | +8.00 | +8.00 |
| +7.00 | +8.00 | +8.00 | +8.00 | +8.00 |

* Table values compensated for vertex distance. No vertex conversion necessary. Shaded boxes represent powers available at a later date.



ALL THE BENEFITS OF A PROGRESSIVE
IN A SOFT CONTACT LENS