

CIBA VISION® Lens Identification Marks

Current Products

AIR OPTIX® sphere	CIBA X6 = 8.6 / 14.2
AIR OPTIX® for ASTIGMATISM	CIBA B7 = 8.7 / 14.5
AIR OPTIX® MULTIFOCAL	CIBA PL = 8.6/14.2,Low add CIBA PM = 8.6/14.2,Med add CIBA PH = 8.6/14.2,High add
AIR OPTIX® NIGHT & DAY® AQUA	CIBA Z4 = 8.4 / 13.8 CIBA Z6 = 8.6 / 13.8
CIBASOFT® PROGRESSIVE TORIC	CIBA CPT86 = 8.6 / 14.5 CIBA CPT89 = 8.9 / 14.5
CIBASOFT® SOFTCOLORS® 13.8	CIBA A3 = 8.3 / 13.8 CIBA A6 = 8.6 / 13.8 CIBA A9 = 8.9 / 13.8
CIBASOFT® SOFTCOLORS 14.5	CIBA B6 = 8.6 / 14.5 CIBA B9 = 8.9 / 14.5 CIBA B2 = 9.2 / 14.5
CIBASOFT® VISITINT® 13.8	CIBA A3 = 8.3 / 13.8 CIBA A6 = 8.6 / 13.8 CIBA A9 = 8.9 / 13.8
CIBASOFT® VISITINT® 14.5	CIBA B6 = 8.6 / 14.5 CIBA B9 = 8.9 / 14.5 CIBA B2 = 9.2 / 14.5
CIBASOFT® STANDARD VISITINT®	CIBA E3 = 8.3 / 13.8 CIBA E6 = 8.6 / 13.8 CIBA E9 = 8.9 / 13.8
CSI® TORIC	C (-1.00 Cylinder) D (-1.75 Cylinder) E (-2.50 Cylinder)
DAILIES® AQUACOMFORT® PLUS	No markings
FOCUS® DAILIES®, FOCUS® DAILIES® PROGRESSIVES, or FOCUS® DAILIES® TORIC	No markings
FOCUS® PROGRESSIVES	CIBA V6 = 8.6 / 14.0 CIBA V9 = 8.9 / 14.0
FOCUS® SOFTCOLORS® MONTHLY	CIBA S6 = 8.6 / 14.0 CIBA S9 = 8.9 / 14.0
FOCUS® VISITINT® MONTHLY	CIBA S6 = 8.6 / 14.0 CIBA S9 = 8.9 / 14.0
FOCUS® TORIC MONTHLY VISITINT®	CIBA M9 = 8.9 / 14.5 CIBA M2 = 9.2 / 14.5
FOCUS® 1-2 VISITINT® (a.k.a. NewVues® prior to '97)	CIBA N4 = 8.4 / 14.0 CIBA N8 = 8.8 / 14.0
FOCUS® 1-2 SOFTCOLORS® (a.k.a. NewVues® prior to '97)	CIBA N4 = 8.4 / 14.0 CIBA N8 = 8.8 / 14.0
FRESHLOOK®	FL
ILLUSIONS®	CIBA I3 = 8.3 / 13.8 CIBA I6 = 8.6 / 13.8 CIBA I9 = 8.9 / 13.8
NIGHT & DAY®	CIBA Z4 = 8.4 / 13.8 CIBA Z6 = 8.6 / 13.8
O₂OPTIX®	CIBA X6 = 8.6 / 14.2
TORISOFT® (-) Powers	CIBA T86 = 8.6 / 14.5 CIBA T89 = 8.9 / 14.5 CIBA T92 = 9.2 / 14.5
TORISOFT® (+) Powers	CIBA +6 = 8.6 / 14.5 CIBA +9 = 8.9 / 14.5 CIBA +2 = 9.2 / 14.5

Discontinued / Sold Products

AOSOFT (Discon. 9/96)	No Markings
AQUAFLEX (Lombart 1/05)	No Markings
BISOFT (Lombart 1/95)	CIBA C3 = 8.3 / 14.0 CIBA C6 = 8.6 / 14.0 CIBA C9 = 8.9 / 14.0
CHOICE A.B.® (Discon. 9/06)	No Markings
CIBASOFT® CLEAR 13.8 (Discon. 1/98)	CIBA A3 = 8.3 / 13.8 CIBA A6 = 8.6 / 13.8 CIBA A9 = 8.9 / 13.8
CIBASOFT® CLEAR 14.5 (Discon. 1/98)	CIBA B6 = 8.6 / 14.5 CIBA B9 = 8.9 / 14.5 CIBA B2 = 9.2 / 14.5
CIBASOFT® STANDARD® CLEAR (Discon. 1/98)	CIBA E3 = 8.3 / 13.8 CIBA E6 = 8.6 / 13.8 CIBA E9 = 8.9 / 13.8
CIBASOFT® STANDARD® SOFTCOLORS® (Discon. 1/98)	CIBA E3 = 8.3 / 13.8 CIBA E6 = 8.6 / 13.8 CIBA E9 = 8.9 / 13.8
CIBATHIN (Discon. 1/98)	CIBA D6 = 8.6 / 13.8 CIBA D9 = 8.9 / 13.8
CIBATHIN SOFTCOLORS® (Discon. 1/98)	D6, D9
PROTEK™ (Discon. 2000)	CIBA P9 or CIBA P2
SOFTCON DW (Lombart 1/95)	No Markings
SOFTCON DW APHAKIC (Discon. 12/92)	No Markings
SOFTCON EW (Lombart 1/05)	No Markings
SOFTCON EW APHAKIC (Lombart 1/95)	No Markings
SPECTRUM BIFOCAL (Discon. 9/96)	CIBA V6 = 8.6 / 14.0 CIBA V9 = 8.9 / 14.0
SPECTRUM SPHERICAL (Discon. 12/94)	CIBA S6 = 8.6 / 14.0 CIBA S9 = 8.9 / 14.0
SPECTRUM TORIC (Discon. 3/95)	CIBA M9 = 8.9 / 14.5 CIBA M2 = 9.2 / 14.5
TORISOFT® SOFTCOLORS® - (Discon. 9/96)	CIBA T86 = 8.6 / 14.5 CIBA T89 = 8.9 / 14.5 CIBA T92 = 9.2 / 14.5
TORISOFT® SOFTCOLORS® + (Discon. 9/96)	CIBA +6 = 8.6 / 14.5 CIBA +9 = 8.9 / 14.5 CIBA +2 = 9.2 / 14.5